

Practical Biotechnology

Lab 5

Competent cell preparation

By Dr. Shaymaa H. Ali

Assistant Professor of Molecular biology

Dr. Bakhtyar N. Ali

PhD in Biotechnology

Mrs. Najat Taher


MSc. In molecular Biology



What are competent cells?

Cell competence refers to a cell's ability to take up foreign (extracellular) DNA from its surrounding environment.

The process of genetic uptake is referred to as *transformation*.



Principle of Competent Cells

Some cells got to be exposed to some **chemical** or **electrical** treatments to transform them into competent cells by altered cell walls that allow the DNA to simply undergo it.


Cell competence has become an essential research tool for cloning because it provides scientist a mechanism to introduce new genetic material into a cell.

Other factors include existing DNA damage within the cell and recombination ability of incoming DNA.




Noncompetent *E. coli*

E. coli that has not been made competent will not take up plasmids introduced into their environment.



Competent *E. coli*

E. coli made competent either through CaCl_2 and heat-shock or through electroporation will have better membrane permeability (pores), enabling plasmid uptake.



Methods of Preparation of Competent Cells

Competence is achieved in two ways:

- Natural Competence
- Artificial Competence

Natural cell competence is genetically determined, that is to say, a bacterium is genetically predisposed to take up free genetic material that exists within their environment.

Artificial or induced competent cells are cells researchers have made competent through electrical (electroporation) or chemical manipulation. **Electroporation** is the process in which cells take up DNA.




Calcium Chloride: The bacterial cells were treated with salt then suddenly exposed to high temperatures. This is known as the heat shock treatment method and it is the standard method for the preparation of those cells.


Electroporation: During this technique, an electrical field is applied to the cells to extend their permeability. It is also known as electroporabilization.



Chemical transformation



Electroporation



Procedure:

1. Inoculate 2ml of LB or (nutrient broth) with a single DH5 α colony. Incubate culture overnight at 37°C while shaking at 250 RPM.

2. The following morning, inoculate 500ml of LB or (nutrient broth) with 1ml of saturated overnight culture. Incubate culture at 37°C while shaking at 250 RPM until OD₆₀₀ = 0.5 (3-5 hours).

3. Transfer culture to 2 pre-chilled sterile 250ml centrifuge tubes.

Pellet bacteria cells with a 5000 RPM spin for 10 minutes at 4°C. Discard supernatant. Place pellets on ice.

4. Resuspend cells in 10ml cold CaCl₂ solution. Pool cells together into one pre-chilled 50ml Oakridge tube.

5. Pellet cells with a 2500 RPM spin for 5 minutes at 4°C. Discard supernatant and resuspend cells in 10ml cold CaCl₂ solution. Set on ice 30 minutes.

6. Pellet cells with a 2500 RPM spin for 5 minutes at 4°C. Discard supernatant and resuspend cells in 2ml cold CaCl₂ solution.

*At this point you can leave cells on ice overnight at 4°C – this increases competency in some cases.

7. Dispense cells into 50ul aliquot in pre-chilled sterile polypropylene tubes.

Store cells at –80°C. Test for Competency Transformation.

Competent Cell Preparation



Resuspend each tube with 40ml chilled CaCl_2 by gentle pipetting

Thank you for
listening



8. Remove competent DH5a cells from the -80°C and immediately place on ice. Once thawed, add $>10\text{ng}$ of plasmid DNA to a $50\mu\text{l}$ aliquot of competent cells. Place cells/DNA on ice for 3 minutes.

9. Heat shock cells at 42°C for 3 minutes.

10. Place cells back on ice for 3 minutes.

11. Add 1ml LB to cells/DNA. Tape tube onto shaking incubator platform and incubate cells/DNA for 1 hour at 37°C while shaking at 250 RPM.

12. Pellet cells with a quick spin. Remove $800\mu\text{l}$ of supernatant. Resuspend cells in the remaining supernatant.

13. Plate $100\mu\text{l}$ and $200\mu\text{l}$ of transformation onto 2 LB+Amp plate. Place plates inverted at 37°C overnight.



Solutions:

1. 100 mM MgCl₂:

- 1:10 dilution of lab stock; use sterile ingredients or filter
- sterilize

2. 100 mM CaCl₂:

- 1:10 dilution of lab stock; use sterile ingredients or filter
- sterilize

3. 85% 100 mM CaCl₂, 15% glycerol:

- 42.5 ml 100 mM CaCl₂
- 7.5 ml 100% glycerol
- 50.0 ml total volume; mix well and use sterile ingredients or filter sterilize