

Practical Medical Bacteriology

Lab 8

Laboratory Diagnosis of Enterobacteriaceae spp: *E. coli* & *Klebsiella*



2021/2022

General Characters of Enterobacteriaceae spp

- Short Gram-negative bacilli (coccobacilli or as straight).
- Their natural habitat is the intestinal tract of humans and animals.
- Many genera like (*Escherichia*, *Shigella*, *Salmonella*, *Enterobacter*, *Klebsiella*, *Serratia*, *Proteus*, and others).
- Most motile with peritrichous flagella except *Shigella*, *Klebsiella* and *Yersinia* are non-motile.
- Non-spore-forming, facultative anaerobic or aerobic microbes.
- Some are opportunistic pathogens and others regularly pathogenic for humans.
- Ferment a wide range of carbohydrates.
- Distinguishing properties associated with all enterobacteriaceae:
Ferment glucose, reduce nitrate to nitrite and oxidase negative



Laboratory diagnosis

1. Specimen:


- *Escherichia coli*: according to disease
 - Urine, pus and stool, pus, blood, stool and CSF
- *Klebsiella*: according to disease
 - Sputum, urine, pus.

2. Microscopic morphology:

- *Escherichia coli*: Gram negative rods, coccobacilli, or straight rods with rounded ends.
- *Klebsiella*: Gram negative rods, with unstained halo around them due to the presence of the capsule.



A
FIGURE 15-1 A: Gram stain of *Escherichia coli*. Original magnification





Laboratory diagnosis

3. Cultural characteristics

❑ On blood agar

- *Escherichia coli*: Colonies are medium sized, smooth, round, grayish-white colonies on blood agar. Some strains produce β hemolysis.
- *Klebsiella*: Slimy appearance of the colonies on blood agar.



Laboratory diagnosis

3. Cultural characteristics

❑ On MacConkey agar

- *Escherichia coli*: deep red colonies are produced as the organism is ferment lactose.
- *Klebsiella*: it shows red/pink mucoid lactose fermenting colonies.





Laboratory diagnosis

3. Cultural characteristics

❑ On Eosin Methylene Blue agar

- *Escherichia coli*: produces black colonies with greenish-black metallic sheen.
- *Klebsiella*: produces large, mucoid, pink to purple colonies with no metallic green sheen on EMB agar.



Laboratory diagnosis

4. Biochemical reactions A. IMViC


Escherichia coli

+ + - -

Indole test

Positive result *E. coli*

Negative result *Klebsiella*




Klebsiella species

- - + +

Methyl Red test

Red: Positive MR (*E. coli*)


Yellow or orange:
Negative MR (*Klebsiella*)



Voges-Proskauer test

Pink: Positive VP
(*Klebsiella*)


No change: Negative VP
(*E. coli*)



Citrate utilization test

Positive *Klebsiella*

Negative *E. coli*




Laboratory diagnosis

4. Biochemical reactions

B. Carbohydrates fermentation in triple sugar iron agar

- *Escherichia coli*: A/A /g+ / H₂S –
- *Klebsiella species*: A/A /g++ / H₂S -



Laboratory diagnosis

4. Biochemical reactions


C. Urease test

- *Escherichia coli*: **Negative**
- *Klebsiella species*: **Positive**



Motility test:

- *Escherichia coli*: **Motile**
- *Klebsiella species*: **Non-motile**



Identification by API20E

